

Paingone & Paingone Plus

1. How does paingone work?

- a) It vibrates and massages the area that it is used on
- b) It emits a warming therapy helping to increase blood flow
- c) It emits a small electrical pulse, which results in the release of endorphins
- d) It emits an infra-red therapy

2. Which one of the following statements regarding Paingone & Paingone Plus is TRUE?

- a) Paingone & Paingone Plus carry a one year warranty against malfunction
- b) Paingone & Paingone Plus carry a two year warranty against malfunction
- c) Paingone & Paingone Plus carry a three year warranty against malfunction
- d) Paingone & Paingone Plus carry a five year warranty against malfunction

3. How many times can you use Paingone Plus during a 24 hour period?

- a) Paingone Plus can be used every 3-4 hours
- b) Paingone Plus can be used up to a maximum of 4 times per day
- c) Paingone Plus can be used as often as required. You cannot overdose
- d) Paingone Plus can be used up to a maximum 10 times per day

4. Which one of the following statements regarding Paingone Plus is TRUE?

- a) Paingone can be used on arthritis, rheumatism, joint pain, tennis elbow & sciatica
- b) Paingone can be used on back pain, osteoporosis, inflammation & migraine
- c) Paingone can be used on most physical pain
- d) all of the above

5. When should you seek medical advice before using paingone?

- a) If you are fitted with a pace maker or suffer from epilepsy
- b) If you are pregnant
- c) If you are taking medication such as steroids/ warfarin
- d) All of the above

6. Which one of the following statements regarding Paingone Plus is TRUE?

- a) paingone Plus is powered by battery
- b) Paingone Plus contains a piezo electric element – no charging is required
- c) Paingone Plus contains a rechargeable battery
- d) Paingone Plus works on main power

7. A customer is interested in the Paingone and Paingone FREEZE cooling gel, what actions would you recommend on using both of them together?

- a) apply the gel before using the Paingone
- b) **apply the gel after using the Paingone**
- c) apply the gel directly before using Paingone
- d) Paingone and Paingone FREEZE cooling gel shouldn't be used together

8. Which one of the following statements regarding Paingone Plus is TRUE?

- a) To use Paingone Plus click to button 30-40 times
- b) To use Paingone Plus press the on/off button, use for 30 seconds then click the on/off button again
- c) **To use Paingone Plus press and hold the grey button for 30 seconds**
- d) None of the above